

Submitting Histology Data

This page provides instructions for adding histological (H&E) images to the GUDMAP/RBK Data Explorer.

If you have any questions or feedback, please send them to your consortium's help email: help@gudmap.org or help@rebuildingakidney.org

We also have the following training materials available:

- [Webinar Slides](#)
- [Webinar Replay \(18:40\)](#)
- Tutorial Videos (Coming Soon)


Overview

Adding histological slides involve the following steps:

- Make sure you are in the correct Globus authentication group, [kidney-writers](#), and that you are logged in.
- Create a new Histological Slide record.
- Fill out the form and upload your CZI file.
- Once the slide is submitted to the Hub (via the *Curation Status* field), the system will process the image.
- The image is then displayed in the browser via a special viewer that displays channels, a scale bar, and allows you to zoom in and out and create annotations.

Schema

Histological Image Schema (Detailed)



The schema for histological slides is very simple - the record for the slide and then the system-generated image.


Are you in the kidney-writers group?

If you haven't already done so, go to this link to join the group: <https://www.globus.org/app/groups/af0b4010-5b75-11e6-9575-22000aef184d/about>

You can find more details about this process at [Accessing GUDMAP and RBK Resources](#).

1. Create Histological Slide record

- In the top navigation bar, click *Create > Histological Images*.

















- The "Create Histological Slide Record" form appears:

RBK/GUDMAP Resources
Search
Dev-Only Search
Create
Feedback
cristina@globusid.org

Create Histological Slides Record
Submit Data
+
v

* indicates required field

Record Number	1		
RID	_ermrest.urlb32_encode(nextval('_ermrest.rid_seq'::regclass))		
* Name	<input type="text"/>		
* Download Link	No file Selected	 Select file	
Tissue	Select a value		
Age Stage	Select a value		
Gender	Select a value		
Specimen Fixation	Select a value		
Embedding Medium	Select a value		
Staining Protocol	Select a value		
* Species	Homo sapiens		
Description	<input type="text"/>		
Pixels Per Meter	<input type="text"/>		
Curation Status	In Preparation		
* Principal Investigator	Select a value		
* Data Provider	Select a value		
* Consortium	Select a value		

- Select the values for each relevant field. The **required** fields are listed, but the more fields you fill out the easier the data is to find:
 - *Name*: Please use a short descriptive name to help users understand what the image represents.
 - *Download Link*: Upload the raw CZI file.
 - *Tissue*
 - *Age Stage*
 - *Specimen Fixation*
 - *Embedding Medium*
 - *Staining Protocol*
 - *Curation Status*: Choose either:
 - *In Preparation*: Use this status while still drafting the data.
 - *PI Review*: Use this status when your data is ready for internal review.
 - *Submitted*: Use this status when your data is ready for Hub review.
 - **Note**: Your data will **not** be viewable publicly until approved for *Release* by the Hub. [For a complete description of the Curation Process, click here.](#)
 - *Principal Investigator*
 - *Data Provider*
 - *Consortium*
- Scroll back up to the top of the page and click the *Submit* button to save the record.

2. Reviewing and Submitting Histological slides

Note: By the hard launch of the new GUDMAP site in April, there will be dashboards and email notifications to make this process more straightforward. In the meantime, here is how a project's PI or designated reviewer can find their project's data with a Curation Status of "PI Review"

- Make sure you are logged in.
- From the navigation bar, click *Search > Gene Expression Data > Histological Images*.
- In the faceting sidebar on the left, scroll to **Curation Status** and choose *PI Review*. Note: Keep in mind that the data submitter may have forgotten to set the Curation Status field, in which case the status would still be *In Preparation*.
- In the faceting sidebar, scroll to **Principal Investigator** and choose your project's PI. Now you should see the data you need to review.
- When your record is approved internally, change *Curation Status* to *Submitted* to send it to the Hub (click here for the full [Curation Workflow](#)).

3. Processing the CZI file

Once you have finished any internal review designated by your lab and you set the *Curation Status* to *Submitted*, the CZI file will be added to the system queue for processing.

The system runs a script every hour that looks in the queue and processes the CZI files there, which will take an additional 30 minutes per file. So it may take up to 1-2 hours for your CZI file to be processed and available in the viewer. The raw file is then immediately available.

Here is an example of a processed CZI file within the viewer:

▼ Images (showing all 1 results)


[HuKid_28_1b \(full screen\)](#)



4. Deleting Histological Slides



To delete the Histological slide record, you must first delete the system-generated image.

- Go to the record you want to delete (see previous section).
- Scroll down to the "Images" section and click "Table Display" (a link to the right side of the section).



- For each image, click the trashcan (Delete) icon in the "Actions" column. A confirmation window will open - click "Confirm".

▼ Images (showing all 1 results)

Actions	Name ↑↓	Light Source Intensity ↑↓	Microscope ↑↓	Mag
 	20160912-hKD6wWT_173-JAM-0-03-000	277 %		

- Once all the images in this table are deleted, scroll back up to the top of the record and click "Delete" in the record header.

